

GENERAL CHEMISTRY LAB COMPONENT CHE101L CONTENT: LAB 3

Dissolution Reactions: Heats of Dissociation

NAME	Mastroots 1	Jofia Artman	SECTION
STUDENT ID	1921079	642	
DATE & TIME	08-10-	-2024	
	22-10-		/ \/ \/ \/ \/ \/ \/ \/ \/ \/ \/ \/ \/ \
NAME OF THE	INSTRUCTOR	Dr. Mohammed	Kobir Uddin (MKUN)
SIGNATURE &	DATE		
REPORT SUBM	IISSION DATE	05-11-2024	

EXPERIMENT 3

SESSION 1 (QUALITATIVE)

DISSOLUTION REACTIONS: HEATS OF DISSOCIATION

Heats (exothermic or endothermic) are associated with chemical reactions. Quantity of heat evolved or absorbed is directly proportional to the amount reacted. Consider the reactions below:

Heat could be generated or absorbed in this reaction. When heat is generated/released from a chemical reaction it is called exothermic reaction (you can feel it by touching the reaction container (warmer) and when heat is absorbed the reaction is called endothermic (colder). When reactions occur in a reaction vessel (e.g., Beaker) in aqueous condition, formation and dissociation of chemical bonds occur simultaneously. Bond formation and dissociation involves heat energy of the system which is expressed by the term Q which is called enthalpy.

PROBLEM STATEMENT: Is heat energy related to chemical reactions, how?

This experiment is subdivided into two parts:

- I. QUALITATIVE
- II. QUANTITATIVE

PART I. QUALITATIVE

DATA COLLECTION:

Place about 30 mL of distilled water into a 50 mL beaker. Suspend a thermometer (having 0.1°C division mark) into the beaker using thermometer clamp and ring stand. Please make sure that the thermometer is not touching the bottom of the beaker, as any movement of the beaker could break the thermometer. Record the temperature of water in the beaker in every 30 seconds for 180 seconds.




FIGURE 1: Experimental setup for dissolution reaction

Place a moderate amount (which would be 1 to 3 cm 3) of supplied anhydrous magnesium sulfate (MgSO $_4$) to the beaker. Mix vigorously with the glass rod for 2 minutes. Record your observations. (2 points)

Repeat this procedure with each of the following compounds with two different amounts (roughly 1:2): (2 points)

- a. Sodium Nitrate, NaNO3 Initial Temp: 26°C, Firal Temp: 19.75°C. So, Endothermic.
- b. Sodium Chloride, NaCl Initial Temp: 26°C, Find Temp: 25.75°C. So, Endo thermic
- c. Hydrated Calcium Chloride, CaCl2.2H2O

 Initial Temp: 26°C, Fired Temp: 29.75°C. So, Enothermic.
- d. Ammonium Nitrate, NH4NO3 Initial Temp: 25°C, Final Temp: 19.75°C. So, Endo therenic E. Mognesium Sulfade, MgSOy

Initial Temp: 26°C, Find Temp: 29.75°C. So, Excothermic

DATA ANALYSIS:

dissolution on treaction

Compare and contrast the behaviors of these compounds. Identify any generalizations that can be made about the chemical reactions observed. What conclusions can be drawn from the data? (4 points)

No NO3, No Clord NH41NO3 enhibit endothermic reactions, meaning they obsorb heat from the surroundings, resulting in a temperature decrease. Co Clo. 2H2O and MySO4 are undergo enothermic reactions, releasing heat into the surroundings and causing a temperature trise. In endothermic reactions Na NO3 and NH4NO3 show significant cooling, while Nach has a minimal deeperature decrease, implying a meation endothermic offed. In enothermic reactions Ca Clo. 2H2O and mySO4 both release heat, causing a temperature increase of enocyty 3.75° in both cases that suggests these compands might release similar amounts of energy. Jodinn-based salts generally exchibit endothermic behaviors. Or the other hand, Calcium and Magnesium salts are enothermic indicating that metal ions are likely to drive reactions that release head. In conclusion, the degree of temperature change gives insight into the strongth of the endo thormic on enothermic reaction. Laptor changes indicate stronger heat obsorpoilion on releasing during

SESSION 2(QUANTITATIVE)

PART II. QUANTITATIVE

DATA COLLECTION:

- a. Accurately weigh a 3 to 5 gm sample of MgSO₄ on the analytical balance. Record the exact mass here. For 4 different trials below measure four different weight samples (e.g., 1,2, 4 & 5 grams respectively).
- b. Suspend the thermometer into a polystyrene cup/coffee cup. Make sure of the thermometer is not touching the bottom of the cup. Measure 20 mL of distilled water by a volumetric cylinder into the cup and stir for 240 second. Record the temperature in every 20 seconds. After 240 seconds add MgSO₄ with vigorous mixing while continuing to record data for 5 minutes.
- c. Determine the temperature change, ΔT , for the reaction. This can be done from the difference of the highest temperature minus the slope of the line go through the points from first 240 seconds of data.
- d. Draw a temperature vs. time graph. Draw the best curve through the points and point out what is happing in each part of the curve.
- e. Repeat the whole procedure with NaNO3.

DATA TABLE 1:

TRAILS

(I) Mass of MgSO ₄ _0.5 gm		(II) Mass of MgSO ₄ _1.0 gm	
Time (s)	Temp(°C)	Time(s)	Temp(°C)
20	25	20	25
40	25	40	25
60	25	60	25
80	26	80	26
100	26	100	26
120	26	120	26
140	26	140	26
160	26.25	160	26.5
180	26.25	180	26.5
200	26.25	200	26.5
220	26.25	220	26.5
240	26.25	240	26.5
260	26.25	260	26.5
280	26.25	280	26.5
300	26.25	300	26.5

(III) Mass of MgSO ₄ _1 <u>.5 gm</u>		(IV) Mass of MgSO _{4_} 2.0 gm	
Time (s)	Temp(°C)	Time(s)	Temp(°C)
20	25.8	20	27
40	25.8	40	27
60	25.8	60	27
80	26	80	28.5
100	2.6	100	28.5
120	26.2	120	28.5
140	26.5	140	28.5
160	26.5	160	28.5
180	26.5	180	29
200	27	200	29
220	27.5	220	29
240	27	240	29
260	26.8	260	29
280	27	280	25
300	27	300	29

△T=27.5°C-25.8°C = 1.7°C. △T=29°C-27°C -2°C.

DATA TABLE 2:

TRAILS

(I) Mass of NaNO ₃ _0.5 gm		(II) Mass of NaNO ₃ _1.00 gm	
Time (s)	Temp(°C)	Time(s)	Temp(°C)
20	25	20	25
40	25	40	25
60	25	60	25
80	24	80	23
100	23.75	100	23
120	23.75	120	22
140	23.75	140	22
160	23.75	160	22
180	24.25	180	22
200	24.50	200	22.
220	24.75	220	22
240	24.75	240	22
260	24.75	260	22
280	24.75	280	22
300	24.75	300	22

(III) Mass of NaNO ₂ _1,5 gm		(IV) Mass of NaNO ₃ _2.0 gm	
Time (s)	Temp(°C)	Time(s)	Temp(°C)
20	26	20	25.5
40	26	40	25.5
60	26	60	25.5
80	26	80	25
100	24	100	25
120	24	120	25
140	24	140	25
160	24	160	24.5
180	23.8	180	24.5
200	23.5	200	24.5
220	23.6	220	24
240	23	240	24
260	23	260	24
280	23	280	24
300	23	300	24

=-3°C

AT= 24°C-25.5°C.

=- 1.5.0

DATA ANALYSIS

 What do you understand form both data sets you recorded and from the other trials? (4 points)

Mysoy is evolvamic, so despondent rises as more 19504 is added indicating heat tralease during dissolution. Nature is arbitharmic, so terporduce drops as Nature distance, absorbing heat from the surroundings. In both compands, increasing the mass leads to a more pronounced thermal effect, but in apposite directions (i.e. trise for Mysoy, full for Nature).

Mysoy shows a stronger enablermic tresponse at higher concentrations, while Nature shows a stronger endohormic affect as the massimereases. The trade of desporative chargest Mysoyis radadicely quick, peaking confly and then stabilizing but the trade of demperature charge of Natures more stadiology but the trade of demperature charge of Natures more stadiology especially with higher masses, indicating that the endothermic proposes that slavely treaching theorems. Myson disolates and traleases heat quickly, reaching theorems equilibrium in less time Natures dissolution absorbs head slavely, with cooling offert stabilizing over a larger porciod. Therefore, both compands atabilize after a period, completion of dissolution and head stansfor



2. Calculate the heat, Q & moles, n, of the reaction both data sets. Take help from the equitation Q = C x M x Δ T. Assume C = 4.18 Joules/gram 0 C and M is the mass of

Mass of moders

Melan mass of 14304= (24+32+ (4×16)) =120g/ml 0.5 N=N03: $\Delta T = 26.25^{\circ}C - 25^{\circ}C = 1.25^{\circ}C$. 0.5 Mysou: 27= 26.25°C - 25°C = 1.25°C. 0=(4.180/gc)x20gx1.25°C=104.50 1 Mg 504; AT = 26.5°C-25°C = 1.5°C 0 = (4.182/g°C) × 20g×1.5°C = 125.40 1.5194504: AT=27.5°C-25.8°C=17.°C. Q= (4.180/g°C) x20gx1.7°C = 142.120

2 MgSay: AT= 29°C-27°C = 2°C.

Q= (4.180/g°c)x20gx 2°C = 167.20

 $\frac{1208/\text{rol}}{1208/\text{rol}} = 0.0042 \text{ mol}$ $\frac{1208/\text{rol}}{1208/\text{rol}} = 0.0083 \text{ mol}$ $\frac{1208/\text{rol}}{1208/\text{rol}} = 0.0125 \text{ mol}$

120g/mol = 0.0167 mel

1Mg NaNOs: AT= 22°C-25°C=-3°C. 0= (4.18]/g°c)×20g×-3°c=-250.8] 1.5NaNoj: AT= 23°C-26°C = - 3°C 0= (4.18]/g°C)×20g×-3°C=-250.8] 2 Mgsoy". AT= 24°C-25.5°C=-1.5°C. 0= (4.180/gc)×20g×-1.5°C=-125.40 :. n1 = 0.59 = 0.0059 mal n=Given Moss : no = 14 = 0.0118 mol n=Girtles Moler Mass : no = 35/mol = 0.0118 mol n=Girtles

:. r3 = 1.59 = 0.0176 mel

:. ny = 28 = 0.0235 mol

3. Plot the collected data as moles, n vs. Q both sets of data. Number of moles can be calculated as n = (mass of sample in gram) / (molecular weight in grams/mole). Try tofind algebraic equations. (4 points)


For Mgsoy, $m = \frac{138 - 125.4}{0.0125 - 0.0083} = 3 \times 10^3$

So, the linear equation is: y=mx+c. = (3×103)x+100

For No Nos, m = -250.9 + 195 0.0176-0.0118 = -9.6×103

So, the linear equation is: y= mn+c= (-9.6×103)×+100

FIGURE: plot here



MENTAL MODEL: Use the chemical equation given above to represent the dissolution reaction in this experiment. Draw a picture(s) which describes what is happening in atomic or in molecular level for either MgSO₄ or NaNO₃ system. How heat release or absorbed can be described from these pictures? (5 points)

Mgsoy:


Chemical equation: Mgsoy(s) water(110) Mg2+ (aq) + SOy- (aq)

When represent sulfat (Mgsoy) is dissolved in hador, it dissociates into Mg2+


(magnesimions) and soy- (sulfate ions).

Hot 000 H+

hoder molecules have hydrogen bonds to hold then dogstor. When the hydrogen bonds breakdown inside nator it requires energy. So it is an endothermic process.



when Mgsoy is added the negative dipole onygen get surveyded to the positive ion and this ordered aroungement of mater releases energy called hydration enthalpy and this process is enothermic process.



The positive dipole of hydrogen surrounds the negative ions and trelease energy. In this case, hydrodionenthology is greater than latice energy. So, the solid structure energy is used up to increase the demporature.

Therefore, the every is released when nator melecules from hydrotion shells arroand the dissociated ions. The attraction between nator melecules and the ions results in the release of every, which is observed as heat in the solution. The overall is enothermic, meaning the

energy treleased from hydrodion of ions enceeds the energy: required to break the ionic bonds. This causes the temporature of the water to trise.

NaNO3:

H (NO3) 1.

chemical equation: Nervoz (5) hador(H2O), Not(or) + NO3-(or)
when sodium nitrade (Nervoz) is dissolved in nator, it dissoites into Net
(Sodium ions) and NO3-(nitrate ions).

Her how have molecules have the hydrogen bonds to hold her hydrogen bonds breekdown inside hat of the hador it requires energy. So, it is on endothermic process.

when NaNo3 is added the negative dipole onygen get sourceunded to the positive ion and that Not I this ordored avaragement of water releases less energy required to break the ionic bonds, resulting in ret energy absorption. This process is called endo thermic process.

The positive dipole of hydrogen survivourds the negative ion and absorred energy. In this case, hydration enthalpy is less than the latice erogy.

Therefore, noder molecules from hydration shells arrand the dissociated ions. However the energy treleased during this hydration process is loss than the energy neguired to brook the ionic bonds, tresulting in net energy obsorphion, making the solution feel cooler. As Nanos disable, morre energy is absorbed than is treleased, leading to a drop in temperature. Thus this reaction is an endo thormic treaction.